

ACTIVE CROWN

Robust in design.
Critical to your production.
Within a trusted ecosystem.



Base

Power	DC 5V @ 3A via USB-C
Mounting	2x 1/4" -20 UNC thread 1x 3/8" -16 UNC thread Compatible with standard-thread mounting plates as listed, including ARRI Pin-Lock and other mounting plates supported by Vicon cameras (mounts not included)
Dimensions	115mm x 72mm
Weight	297g (including complete set of IR LED Stalks)

IR LED Stalks

Quantity	6x (3x small, 1x medium, 2x large)
Markers	Active LED 150° Field of View Always on (no sync required)
Dimensions	Small - 30mm x 27mm x 80mm Medium - 30mm x 27mm x 159mm Large - 30mm x 27mm x 179mm
Mounting	Magnetic sockets on base
Visibility	VK26 can see the crown at 50m (*metrics on all release cameras are in review)

Peripherals

Power	1 x 1m USB-C to USB-C cable
Other	7x blanking caps for unused stalk sockets

www.vicon.com/activecrown

VICON